# Instructions For Use

## Blue Sky Bio, LLC Dental Implant Systems

<b>English</b>	Blue Sky Bio, LLC Dental Implant Systems
<b>Le français</b>	Systèmes d'implants dentaires Blue Sky Bio, LLC
<b>Deutsch</b>	Blue Sky Bio, LLC-Dentalimplantatsysteme
<b>Italiano</b>	Sistemi di impianto dentale Blue Sky Bio, LLC
<b>Español</b>	Sistemas de implantes dentales Blue Sky Bio, LLC
<b>Svenska</b>	Blue Sky Bio, LLC Dental Implant Systems
<b>Türkçe</b>	Blue Sky Bio, LLC Dental Implant Sistemleri
<b>Ελληνικά</b>	Συστήματα οδοντικών εμφυτευμάτων Blue Sky Bio, LLC
<b>Português</b>	Sistemas de implante dentário da Blue Sky Bio, LLC

**Blue Sky Bio, LLC**  
800 Liberty Dr  
Libertyville, IL 60048  
USA

EC REP

**mdi Europa GmbH**  
Langenhagener Straße 71  
D-30855 Langenhagen  
Germany

CE 0297

REF

CATALOG NUMBER

DO NOT REUSE

CAUTION, CONSULT ACCOMPANYING DOCUMENTS

Rx Only

CAUTION FEDERAL LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN OR PROPERLY LICENSED PRACTITIONER

MANUFACTURED BY

STERILE R

STERILIZED USING GAMMA RADIATION

DO NOT USE IF PACKAGE IS DAMAGED

LOT

BATCH CODE

EC REP

AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY

USE BY DATE

**elfu Indicator**

ELECTRONIC INSTRUCTIONS FOR USE

NON-STERILE

MD

MEDICAL DEVICE

UDI

UNIQUE DEVICE IDENTIFIER

DATE OF MANUFACTURE

SINGLE STERILE BARRIER SYSTEM WITH PROTECTIVE PACKAGING INSIDE

# English (en) (For US Only)

## Blue Sky Bio, LLC Dental Implant Systems

### Instructions For Use

#### Indications

The implant systems are intended for use in supporting single or multiple restorations in the fully or partially edentulous mandibular or maxillary alveolar process.

#### Contraindications

Implants should not be placed anytime when there are general contraindications associated with elective oral surgery. Absolute and relative contraindications include, but are not limited to: Titanium allergy (Chrome, Cobalt or Molybdenum allergy for abutments manufactured from such alloys), cardiac and vascular disease, bleeding disorders, psychological disorders, uncontrolled diabetes mellitus, mineral, bone, or connective tissue disorders, renal disease, hepatic disease, auto-immune disorders, decreased immune function due to disease or medications, infectious disorders, and adverse conditions caused by medications. Further relative contraindications include poor oral hygiene, bruxism, malnutrition, alcoholism, tobacco usage, and history of radiation therapy. In addition, the patient needs an adequate volume of residual bone for the placement of implants of sufficient size and number to support the anticipated functional loads to which the patient will subject these implants.

Narrow implants and angled abutments are not intended for use in the posterior region of the mouth.

#### Warning

Implants should be placed and restored only by practitioners who are licensed and trained to perform these procedures. Adequate preoperative studies should be performed to examine the anatomic structures and to assess the biomechanical, functional, and esthetic requirements of each case. Radiographs or other diagnostic reviews should be performed to determine position and topography of the maxillary sinus, nasal cavities, inferior alveolar nerve, mental foramen, natural tooth positions and other anatomical features that may affect implant placement or prognosis. Consultation between the surgeon, restorative dentist, and dental laboratory is essential for success. Risks of implant placement and restoration include, but are not limited to: infection, implant failure, loss of bone and soft tissue, unfavorable aesthetic result, anesthesia, dysesthesia and paresthesia in the oral and facial areas, sinus infection, dislodgement of implants and instruments in the surrounding structures, damage to adjacent teeth, non-restorable implants, fracture of implants or restorative components, and loosening of implants or restorative components. Each implant system has unique measuring characteristics to allow full seating of the implant to the desired depth. In some instances, drill length reference lines measure longer than the stated length of the implant. It is recommended that the implant surgeon be thoroughly familiar with the specific measurement system being utilized and provide a suitable safety margin adjacent to any teeth and vital structures. Failure to recognize the difference between the actual length of the drill and radiographic measurements can result in permanent injury to the nerves or other vital structures by drilling beyond the depth intended, potentially resulting in permanent numbness to the lower lip and chin or other injuries. Each implant system has specific design characteristics for mating implants, abutments, prosthetic components, and instrumentation. Combining instruments and components that are not configured or dimensioned for correct mating can lead to mechanical failure of components, damage to tissue, or unsatisfactory aesthetic results therefore one-hundred percent success cannot be guaranteed. Lack of adequate quantity and/or quality of remaining bone, infection, inadequate surgical technique, poor patient oral hygiene, and generalized disease are some potential causes for failure of osseointegration, both immediately after surgery or after osseointegration is initially achieved. Pre-operative hard tissue or soft tissue deficits may yield a compromised aesthetic result or unfavorable implant angulation. With respect to children, treatment is not recommended until completion of alveolar growth has been verified.

#### MRI Safety Information

While these products have not undergone MRI safety testing, literature analysis indicates that scanning a Blue Sky Bio, LLC Implant System poses no significant risks under specific conditions and are considered as "MR Conditionally Safe" under the following scanning conditions:

- Static magnetic field strengths of  $\leq 3$  Tesla
- Magnetic field gradients up to 30T/M (3000G/cm).
- RF Excitation:  
Circularly Polarized (CP) for all landmarks above umbilicus - No excitation restrictions for all landmarks below the umbilicus
- Operating Mode: Normal Operating Mode in the allowed imaging zone
- Whole-body specific absorption rates (SAR) of 2W/kg
- Maximum Head SAR: Not evaluated for head landmark
- Scan Duration: No specific constraints due to implant heating

#### Procedural Precautions, Surgery

All efforts must be made to minimize damage to the host tissue. In particular, special attention must be paid to thermal and surgical trauma and to the elimination of contaminants and sources of infection. The surgical procedure requires a high degree of precision and care. Any divergence from the principle of least possible trauma at implant installation increases the risk of failure to establish osseointegration. All drilling procedures should be performed at maximum 1000-2000 RPM with copious irrigation. The use of sharp drills, sufficient irrigation, an in-and-out drilling motion, short cutting cycles, waiting for the bone to cool, and use of pilot drills in successively increasing sizes are essential. Please refer to our web site for the specific sequence of drills for each implant type and size. Special care should be taken to avoid over or under preparation of the osteotomy. Implants should be inserted in such a way that they are stable and lack any mobility. Excessive insertion torque (greater than 60 Ncm) may lead to damage to the implant or instruments and fracture or necrosis of the bone site. All instruments used in surgery must be maintained in good condition and care must be taken that the instruments do not damage the implants or other components. Precautions must be taken to avoid the swallowing or aspiration of components used in implant dentistry. After the implant installation, the surgeon's evaluation of bone quality and initial stability will determine when implants may be loaded. An appropriate follow-up protocol should be followed.

#### Procedural Precautions, Prosthetics

Especially important is proper stress distribution through passive adaptation and fitting of the bridge to the implant abutments, adjusting occlusion to the opposing jaw, and avoiding excessive transverse loading forces. Immediate loading and immediate temporization require additional precautions and are not suitable for all cases. Because of the small size of prosthetic components, care must be taken that they are not swallowed or aspirated by the patient.

#### Instructions For Abutments For Restoration Of Blue Sky Bio, LLC Implants

After adequate osseointegration the implants should be uncovered, if necessary, and the cover screw should be removed. The internal part of the implant should be irrigated, freed from debris, and dried. If necessary, a healing abutment should be placed and the tissue should be allowed to heal around the healing abutment. A final abutment that is compatible with the implant should be chosen. If the abutment is modified at chair side, it is advisable to make an impression and make a preliminary model to identify any undercuts prior to making the final impression. A carbide or diamond bur can be used with copious irrigation to remove the undercuts. The indexing feature of the abutment (e.g. octagon or trilobe) should match the internal indexing feature of the implant. Adequate seating can be verified with a radiograph.

If the abutment is modified in a laboratory, proper orientation of the abutment into the implant should be identified by means of a transfer jig or any other transfer device. Final abutment seating should be performed using an insertion driver compatible with the selected implant system and a torque wrench applying 30 Ncm of torque on the abutment or fastening screw. Abutments should not be over-prepared, and the retention of restoration and the strength of the abutment should be taken into consideration. It is recommended that less than 50% of the structure of customizable abutments and less than 30% of straight abutments be removed. Do not modify abutment at the junction with the implant. The abutment angle should not be changed from the specified angulation. Substantial modification of the abutment may compromise the structural integrity and can lead to premature failure. Do not modify abutment at the junction with the implant. A restoration should be fabricated and inserted, while making sure that the restoration is stable and occlusal load is appropriate. Multi-unit abutments for fixed restorations are to be used only for multiple splinted restorations.

## **Instructions for Base Abutments and for Custom abutments from Titanium Preformed Blanks**

### **Indications for use**

Blue Sky Bio CAD-CAM Abutments are intended to be used in conjunction with Blue Sky Bio endosseous dental implants in the maxillary or mandibular arch to provide support for single-unit or multi-unit prosthetic restorations.

All digitally designed abutments for use with Blue Sky Bio CAD-CAM Abutments are intended to be sent to a Blue Sky Bio validated milling center for manufacture.

### **Description**

Blue Sky Bio CAD-CAM Abutments include Titanium Base Abutments, Cobalt Base Abutments, and Titanium Blank Abutments. Titanium Base Abutments and Cobalt Base Abutments are provided in straight and 15° angled designs. Titanium Base Abutments and Cobalt Base abutments can support a zirconia superstructure or a final restoration. Cobalt Base Abutments may be cast to a final one-piece abutment.

### **Fabrication of Superstructure for Titanium and Cobalt Base Abutments**

The Blue Sky Bio validated milling center is sent the intraoral scan STL file and designs the superstructure using a computer software (CAD) the according to the design parameters below and the interproximal space available.

The recommend material for milling (CAM) the superstructure is Zirconia Y-TZP, conforming to ISO 13356 Implants for surgery - Ceramic materials based on yttria-stabilized tetragonal zirconia (Y-TZP).

### **Superstructure Bonding**

The superstructure may be bonded to the abutment by the milling center, or may be bonded by the clinician. The following instructions are to be followed in either case.

Inspect the abutment for any defects. If any visible defects are present, do not use the abutment.

Grit blasting of the base abutment is not required for proper bonding due to the retentive features on the base.

Cleaning the base and the superstructure. Confirm there is no debris, grease, or other foreign matter present on any area to be cemented.

Multilink Hybrid Abutment Cement (Ivoclar Vivadent AG) is recommended. Refer to the manufacturer's instructions for proper mixing.

Apply cement to the contact surfaces of the base abutment and press the superstructure onto the abutment following the anti-rotational features.

Verify that both components are seated completely.

Remove excess cement immediately and refer to the manufacturer's instructions for proper curing procedures.

Sterilize abutment according to the steam autoclaving instructions.

### **Cobalt Base Abutments Workflow for Casting**

**Cobalt Base Abutments (Straight) can be used as a base for a cast-to straight abutments only.**

**Angled Cobalt Base Abutments are not indicated for cast-to final abutment).**

The final one-piece abutment is fabricated using standard lost wax casting technique.

### **Blue Sky Bio Titanium Preformed Blank Abutments**

The Blue Sky Bio validated milling center is sent the intraoral scan STL file and designs the final abutment using a computer software (CAD) the according to the design parameters below and the interproximal space available. The abutment is machined (CAM) using standard milling techniques.

#### **Design Parameters**

The design parameters for Titanium Base Abutments, Cobalt Base Abutments (including cast-to abutments) and for custom abutments fabricated from Titanium Preformed Blank Abutments are:

Minimum wall thickness – 0.4 mm

Minimum post height for single-unit restorations – 4.0 mm;

Minimum gingival height – 0.5 mm

Maximum gingival height – 6.7 mm

#### **Maximum angulation**

Base Abutments (Straight) – superstructure maximum angle – 30°

Base Abutments (Angled 15°) – superstructure maximum angle – 15° (total of 30°)

Final custom milled abutments – 30°

Cobalt cast-to abutments - 0°

Maximum divergence or angular correction between two implants is 60°

### **Warning**

It is contraindicated to fabricate a zirconia superstructure and to cast a one-piece abutment, or to mill an abutment that exceed the design parameters.

### **Follow Up Care**

Patients should be instructed in appropriate oral hygiene and care of the implants and restorations. Periodic follow up appointments should be made to confirm and maintain adequate

function of the implants and the health of the surrounding tissues.

**Single Use Products:**

Single Use products must not be reprocessed. Reprocessing could cause loss of mechanical, chemical and/or biological characteristics. Reuse could cause cross contamination.

Single-use instruments are marked on the label by the "Do not reuse" symbol and should be disposed of after surgical usage.

Single-use instruments must not be reprocessed after surgical use. Re-use of single-use instruments creates a potential risk of contamination or device failure, which may lead to injury, illness or death of the patient..

**Reusable Instruments:**

Frequent reprocessing has minor effects on the instruments. The end of the product life is determined by wear and damage during use. Therefore, general instruments can be reused with appropriate care, provided they are undamaged and not contaminated

All parts of the instruments must be carefully assessed for visible damage and/or proper functioning.

Device should be visually inspected for any damage prior to each use. Do not use the device if any damage is observed

**Sterility**

All implants are supplied sterile, and are for single use only prior to the labeled expiration date (if applicable). Do not use implants if the packaging has been damaged or previously opened. Abutments and instruments are supplied non-sterile and must be cleaned and sterilized prior to use.

Medical devices that are supplied non-sterile must be sterilized as follows before they are used in the patient:

**Steam Sterilization procedure:** Place the autoclave pouch containing item to be sterilized into the autoclave and follow the specific instructions provided by the manufacturer for pouched items. In general, a pouch must be sterilized by heating for 30 minutes at 250 deg F (121 deg C).

**Cleaning instructions for reusable instruments**

Rinse with cool-to-lukewarm water for two-and-one-half minutes. For all parts place in an ultrasonic cleaner with an enzymatic detergent diluted with tap water per the manufacture's guidelines. Sonicate for 10 minutes. Rinse with tap water for three minutes.

**Method of Supply**

Blue Sky Bio, LLC implants are made out of medical grade titanium alloy (ASTM F136).

Blue Sky Bio, LLC abutments are made out of medical grade titanium alloy , ASTM F136 , chrome-cobalt ASTM F1537 alloy or gold alloy.

**Caution**

The sale of this device is restricted to, or by the order of, licensed physicians or dentists.

**Compatibility:**

BIO   Max NP	BIO   Max NP is surgically and prosthetically compatible with the NobelActive NP implant system.
BIO   Max RP	BIO   Max RP is surgically and prosthetically compatible with the NobelActive RP implant system.
BIO   Internal Hex	BIO   Internal Hex is surgically and prosthetically compatible with the Zimmer/ZimVie Tapered Screw Vent and Implant Direct Internal Hex Legacy3 implant systems.
BIO   Quattro	BIO   Quattro is surgically and prosthetically compatible with the Straumann Bone Level implant system.
BIO   One Stage	BIO   One Stage is surgically and prosthetically compatible with the Straumann Tissue Level implant system.
BIO   Conus	BIO   Conus is surgically and prosthetically compatible with the Astra OsseoSpeed TX implant system.

**Disposal**

The product should be disposed of in an environmentally friendly manner in accordance with local regulations. Hazardous waste of contaminated devices or sharp objects should be disposed of in suitable containers which meet the specific technical requirements.

**Implant Card Instructions for Completion**

An implant card has been provided with the implant device. It is intended to be completed by the healthcare provider; then giving to the patient.

1. Name of the patient or patient ID – To be completed by healthcare institution/provider
2. Date of implantation - To be completed by healthcare institution/provider
3. Name and address of the healthcare institution and provider - To be completed by healthcare institution/provider
4. Barcode to patient implant information – Present on Implant Card

**BlueSkyBio** International Implant Card

 \_\_\_\_\_ **1**

 \_\_\_\_\_ **2**

 \_\_\_\_\_ **3**

\_\_\_\_\_

\_\_\_\_\_

 [www.blueskybio.com/pages/patientimplantinfo](http://www.blueskybio.com/pages/patientimplantinfo) **4** 

**BlueSkyBio** International Implant Card

 \_\_\_\_\_ **Jane Doe**

 \_\_\_\_\_ **03/15/2021**

 \_\_\_\_\_ **XYZ Surgical Center**

\_\_\_\_\_ **123 Healthcare Lane, Chicago IL**

\_\_\_\_\_ **Dr. A.B. Curtis**

 [www.blueskybio.com/pages/patientimplantinfo](http://www.blueskybio.com/pages/patientimplantinfo) 

Locations of text to be filled on the Left. Example of a completed implant card on the right with written information in blue.

Symboling located on the device specific side of the implant card can be found on the first page of the IFU.

Please report any serious incident to the manufacturer

 Blue Sky Bio, LLC  
 800 Liberty Drive  
 Libertyville, IL 60048  
 USA

# English (en) (For UK, Europe, Australia, and Canada)

## Blue Sky Bio, LLC Dental Implant Systems

### Instructions For Use

#### Indications

The implant systems are intended for use in supporting single or multiple restorations in the fully or partially edentulous mandibular or maxillary alveolar process.

#### Contraindications

Implants should not be placed anytime when there are general contraindications associated with elective oral surgery. Absolute and relative contraindications include, but are not limited to: Titanium allergy, cardiac and vascular disease, bleeding disorders, psychological disorders, uncontrolled diabetes mellitus, mineral, bone, or connective tissue disorders, renal disease, hepatic disease, auto-immune disorders, decreased immune function due to disease or medications, infectious disorders, and adverse conditions caused by medications. Further relative contraindications include poor oral hygiene, bruxism, malnutrition, alcoholism, tobacco usage, and history of radiation therapy. In addition, the patient needs an adequate volume of residual bone for the placement of implants of sufficient size and number to support the anticipated functional loads to which the patient will subject these implants.

Narrow implants and angled abutments are not intended for use in the posterior region of the mouth.

#### Warning

Implants should be placed and restored only by practitioners who are licensed and trained to perform these procedures. Adequate preoperative studies should be performed to examine the anatomic structures and to assess the biomechanical, functional, and esthetic requirements of each case. Radiographs or other diagnostic reviews should be performed to determine position and topography of the maxillary sinus, nasal cavities, inferior alveolar nerve, mental foramen, natural tooth positions and other anatomical features that may affect implant placement or prognosis. Consultation between the surgeon, restorative dentist, and dental laboratory is essential for success. Risks of implant placement and restoration include, but are not limited to: infection, implant failure, loss of bone and soft tissue, unfavorable aesthetic result, anesthesia, dysesthesia and paresthesia in the oral and facial areas, sinus infection, dislodgement of implants and instruments in the surrounding structures, damage to adjacent teeth, non-restorable implants, fracture of implants or restorative components, and loosening of implants or restorative components. Each implant system has unique measuring characteristics to allow full seating of the implant to the desired depth. In some instances, drill length reference lines measure longer than the stated length of the implant. It is recommended that the implant surgeon be thoroughly familiar with the specific measurement system being utilized and provide a suitable safety margin adjacent to any teeth and vital structures. Failure to recognize the difference between the actual length of the drill and radiographic measurements can result in permanent injury to the nerves or other vital structures by drilling beyond the depth intended, potentially resulting in permanent numbness to the lower lip and chin or other injuries. Each implant system has specific design characteristics for mating implants, abutments, prosthetic components, and instrumentation. Combining instruments and components that are not configured or dimensioned for correct mating can lead to mechanical failure of components, damage to tissue, or unsatisfactory aesthetic results therefore one-hundred percent success cannot be guaranteed. Lack of adequate quantity and/or quality of remaining bone, infection, inadequate surgical technique, poor patient oral hygiene, and generalized disease are some potential causes for failure of osseointegration, both immediately after surgery or after osseointegration is initially achieved. Pre-operative hard tissue or soft tissue deficits may yield a compromised aesthetic result or unfavorable implant angulation. With respect to children, treatment is not recommended until completion of alveolar growth has been verified.

#### MRI Safety Information

While these products have not undergone MRI safety testing, literature analysis indicates that scanning a Blue Sky Bio, LLC Implant System poses no significant risks under specific conditions and are considered as "MR Conditionally Safe" under the following scanning conditions:

- Static magnetic field strengths of  $\leq 3$  Tesla
- Magnetic field gradients up to 30T/M (3000G/cm).
- RF Excitation:  
Circularly Polarized (CP) for all landmarks above umbilicus - No excitation restrictions for all landmarks below the umbilicus
- Operating Mode: Normal Operating Mode in the allowed imaging zone
- Whole-body specific absorption rates (SAR) of 2W/kg
- Maximum Head SAR: Not evaluated for head landmark
- Scan Duration: No specific constraints due to implant heating

#### Procedural Precautions, Surgery

All efforts must be made to minimize damage to the host tissue. In particular, special attention must be paid to thermal and surgical trauma and to the elimination of contaminants and sources of infection. The surgical procedure requires a high degree of precision and care. Any divergence from the principle of least possible trauma at implant installation increases the risk of failure to establish osseointegration. All drilling procedures should be performed at maximum 1000-2000 RPM with copious irrigation. The use of sharp drills, sufficient irrigation, an in-and-out drilling motion, short cutting cycles, waiting for the bone to cool, and use of pilot drills in successively increasing sizes are essential. Please refer to our web site for the specific sequence of drills for each implant type and size. Special care should be taken to avoid over or under preparation of the osteotomy. Implants should be inserted in such a way that they are stable and lack any mobility. Excessive insertion torque (greater than 60 Ncm) may lead to damage to the implant or instruments and fracture or necrosis of the bone site. All instruments used in surgery must be maintained in good condition and care must be taken that the instruments do not damage the implants or other components. Precautions must be taken to avoid the swallowing or aspiration of components used in implant dentistry. After the implant installation, the surgeon's evaluation of bone quality and initial stability will determine when implants may be loaded. An appropriate follow-up protocol should be followed.

#### Procedural Precautions, Prosthetics

Especially important is proper stress distribution through passive adaptation and fitting of the bridge to the implant abutments, adjusting occlusion to the opposing jaw, and avoiding excessive transverse loading forces. Immediate loading and immediate temporization require additional precautions and are not suitable for all cases. Because of the small size of prosthetic components, care must be taken that they are not swallowed or aspirated by the patient.

#### Instructions For Abutments For Restoration Of Blue Sky Bio, LLC Implants

After adequate osseointegration the implants should be uncovered, if necessary, and the cover screw should be removed. The internal part of the implant should be irrigated, freed from debris, and dried. If necessary, a healing abutment should be placed and the tissue should be allowed to heal around the healing abutment. A final abutment that is compatible with the implant should be chosen. If the abutment is modified at chair side, it is advisable to make an impression and make a preliminary model to identify any undercuts prior to making the final impression. A carbide or diamond bur can be used with copious irrigation to remove the undercuts. The indexing feature of the abutment (e.g. octagon or trilobe) should match the internal indexing feature of the implant. Adequate seating can be verified with a radiograph.

If the abutment is modified in a laboratory, proper orientation of the abutment into the implant should be identified by means of a transfer jig or any other transfer device. Final abutment seating should be performed using an insertion driver compatible with the selected implant system and a torque wrench applying 30 Ncm of torque on the abutment or fastening screw. Abutments should not be over-prepared, and the retention of restoration and the strength of the abutment should be taken into consideration. It is recommended that less than 50% of the structure of customizable abutments and less than 30% of straight abutments be removed. Do not modify abutment at the junction with the implant. The abutment angle should not be changed from the specified angulation. Substantial modification of the abutment may compromise the structural integrity and can lead to premature failure. Do not modify abutment at the junction with the implant. A restoration should be fabricated and inserted, while making sure that the restoration is stable and occlusal load is appropriate. Multi-unit abutments for fixed restorations are to be used only for multiple splinted restorations.

**Follow Up Care**

Patients should be instructed in appropriate oral hygiene and care of the implants and restorations. Periodic follow up appointments should be made to confirm and maintain adequate function of the implants and the health of the surrounding tissues.

**Single Use Products:**

Single Use products must not be reprocessed. Reprocessing could cause loss of mechanical, chemical and/or biological characteristics. Reuse could cause cross contamination.

Single-use instruments are marked on the label by the "Do not reuse" symbol and should be disposed of after surgical usage.

Single-use instruments must not be reprocessed after surgical use. Re-use of single-use instruments creates a potential risk of contamination or device failure, which may lead to injury, illness or death of the patient..

**Reusable Instruments:**

Frequent reprocessing has minor effects on the instruments. The end of the product life is determined by wear and damage during use. Therefore, general instruments can be reused with appropriate care, provided they are undamaged and not contaminated

All parts of the instruments must be carefully assessed for visible damage and/or proper functioning.

Device should be visually inspected for any damage prior to each use. Do not use the device if any damage is observed.

**Sterility**

All implants are supplied sterile, and are for single use only prior to the labeled expiration date (if applicable). Do not use implants if the packaging has been damaged or previously opened.

Abutments and instruments are supplied non-sterile and must be cleaned and sterilized prior to use.

Medical devices that are supplied non-sterile must be sterilized as follows before they are used in the patient:

**Steam Sterilization procedure:** Place the autoclave pouch containing item to be sterilized into the autoclave and follow the specific instructions provided by the manufacturer for pouched items. In general, a pouch must be sterilized by heating for 30 minutes at 250 deg F (121 deg C).

**Cleaning instructions for reusable instruments**

Rinse with cool-to-lukewarm water for two-and-one-half minutes. For all parts place in an ultrasonic cleaner with an enzymatic detergent diluted with tap water per the manufacture's guidelines. Sonicate for 10 minutes. Rinse with tap water for three minutes.

**Method of Supply**

Blue Sky Bio, LLC implants and abutments are made out of medical grade titanium (ASTM F136) or gold alloy.

**Caution**

The sale of this device is restricted to, or by the order of, licensed physicians or dentists.

**Compatibility:**

BIO   Max NP	BIO   Max NP is surgically and prosthetically compatible with the NobelActive NP implant system.
BIO   Max RP	BIO   Max RP is surgically and prosthetically compatible with the NobelActive RP implant system.
BIO   Internal Hex	BIO   Internal Hex is surgically and prosthetically compatible with the Zimmer/ZimVie Tapered Screw Vent and Implant Direct Internal Hex Legacy3 implant systems.
BIO   Quattro	BIO   Quattro is surgically and prosthetically compatible with the Straumann Bone Level implant system.
BIO   One Stage	BIO   One Stage is surgically and prosthetically compatible with the Straumann Tissue Level implant system.
BIO   Conus	BIO   Conus is surgically and prosthetically compatible with the Astra OsseoSpeed TX implant system.

**Disposal:**

The product should be disposed of in an environmentally friendly manner in accordance with local regulations. Hazardous waste of contaminated devices or sharp objects should be disposed of in suitable containers which meet the specific technical requirements.

**Implant Card Instructions for Completion**

An implant card has been provided with the implant device. It is intended to be completed by the healthcare provider; then giving to the patient.

1. Name of the patient or patient ID – To be completed by healthcare institution/provider
2. Date of implantation - To be completed by healthcare institution/provider
3. Name and address of the healthcare institution and provider - To be completed by healthcare institution/provider
4. Barcode to patient implant information – Present on Implant Card

**BlueSkyBio** International Implant Card

 \_\_\_\_\_ 1

 \_\_\_\_\_ 2

 \_\_\_\_\_ 3

\_\_\_\_\_

\_\_\_\_\_

 [www.blueskybio.com/pages/patientimplantinfo](http://www.blueskybio.com/pages/patientimplantinfo)  4

**BlueSkyBio** International Implant Card

 Jane Doe

 03/15/2021

 XYZ Surgical Center

123 Healthcare Lane, Chicago IL

Dr. A.B. Curtis

 [www.blueskybio.com/pages/patientimplantinfo](http://www.blueskybio.com/pages/patientimplantinfo) 

Locations of text to be filled on the Left. Example of a completed implant card on the right with written information in blue.

Symboling located on the device specific side of the implant card can be found on the first page of the IFU.

**For customers in the European Union:**

To request a copy of the Summary of Safety and Clinical Performance (SSCP) for implantable devices, please send an email specifying the article number(s) to [info@blueskybio.com](mailto:info@blueskybio.com)

Please report any serious incident to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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