

















# I-MAX 3D

## 2D programs

<p>PANORAMIC Programs range</p> 	<ul style="list-style-type: none"> <li> L Standard Panoramic</li> <li> L Low dose Panoramic</li> <li> L Half Panoramic Left</li> <li> L Half Panoramic Right</li> <li> L Frontal dentition Panoramic</li> <li> L Improve dentition Panoramic</li> </ul>
<p>BITEWING Programs range</p> 	<ul style="list-style-type: none"> <li> L Bilateral Bitewing</li> <li> L Bitewing Left</li> <li> L Bitewing Right</li> </ul>
<p>SINUS program</p> 	<p>Sinus Standard</p>
<p>TMJ (Temporal Mandibular Joint) Programs range</p> 	<ul style="list-style-type: none"> <li> x4 Standard TMJ</li> <li> x2 Single TMJ</li> </ul>
<p>TEST MODE</p> 	<div style="border: 1px solid black; padding: 5px; display: inline-block;"> <p style="font-size: 2em; margin: 0;">2D</p> </div> <p>Test mode for 2D</p>

**3D programs**

3D Programs range



Full volume: 12x10 cm








Full volume: 9x9 cm




Half volume : 9x5 cm




Small volume: 5x5 cm



Scan of 3D models



radiological guides



impression trays



plasters models

TEST MODE

TEST	<input checked="" type="checkbox"/>	_____
	<input checked="" type="checkbox"/>	_____
	<input checked="" type="checkbox"/>	_____
	<input checked="" type="checkbox"/>	_____



Test mode for 3D